



CAPITAL MEDICAL EQUIPMENTS SALES & SERVICES

MediworldIndia

contact@mediworldindia.com
mediworld21@gmail.com
+91 98900 67867 or +91 8668334988

Bariatric Operating Table

DESCRIPTION

This Electromatic Operating Table especially designed with centric column to bear maximum patient weight capacity of 300 kgs. Its a mobile, electrical operated surgical table for multipurpose i.e. cardiology, bariatric, neurology, laparoscopy & OB/Gynae with the addition of table accessories. The table is equipped with a large translucent sliding table top providing Laparoscopy & C-ARM functions without patient reversing. The Surgical Table is powered by either internal battery or facility electric power system

Features

- Specialized for bariatric, cardiology, OB/Gynae, neurology and laparoscopy.
- Offering full imaging facility from heart to toes by means of C-Arm access.
- With our tabletop longitudinal slide feature, there is no need to reverse patient to achieve full-body fluoroscopic procedures.
- Smooth functions of the table (height, trendelenburg, lateral tilt, flex-reflex, chair, longitudinal top) can be achieved by using the corded handset.
- Table is made of high quality medical stainless steel 304G. Non-reflecting surface, antibacterial, anti-rusted & is easy to clean.
- Interchangeable head/ leg sections allow unobstructed.
- Inbuilt battery enables to run table in emergency



CAPITAL MEDICAL EQUIPMENTS SALES & SERVICES

MediworldIndia

contact@mediworldindia.com
mediworld21@gmail.com
+91 98900 67867 or +91 8668334988



MediworldIndia

contact@mediworldindia.com

mediworld21@gmail.com

+91 98900 67867 or +91 8668334988

CAPITAL MEDICAL EQUIPMENTS SALES & SERVICES

Optional Features

- **Override Panel/Standby Controller:** Additional control panel on the table column cover for continue operation in case of remote failure.
- **Zero Position** to make table neutral by simply pressing '0' positioning button.
- **Memory Function** – The microprocessor based control system is fully programmable and can store upto **two user-programmable preset positions** in its memory that can be recalled anytime by simply pressing M1 or M2 button.
- Remote has **LED backlit** with symbolic position figures to use even in darkness.
- **Safety key** to block motor system in critical operations
- Individual Locking Function

Technical Specifications

Top dimension	L 1960 x W 533-633 mm
Height adjustment	650 mm – 800 mm
Table top sliding	250 mm
Trendelenburg / Reverse	30° / 25°
Lateral tilt	20° / 20°
Kidney elevator	150 mm
Back Rest (up / down)	80° / 25°
Leg Rest (up / down)	15° / 90°
Head Rest (up / down)	20° / 60°
Power supply	24VDC
Battery Backup	3 hours
Patient Weight Capacity	300 kgs (666 lbs)



MediworldIndia

contact@mediworldindia.com

mediworld21@gmail.com

+91 98900 67867 or +91 8668334988

CAPITAL MEDICAL EQUIPMENTS SALES & SERVICES