



CAPITAL MEDICAL EQUIPMENTS SALES & SERVICES

**MediworldIndia**

contact@mediworldindia.com  
mediworld21@gmail.com  
+91 98900 67867 or +91 8668334988

## C-ARM Gearmatic Hydraulic Operation Table

### Description

This is hydraulic surgical table designed especially for Orthopedic, Neurology, Urology and ERCP Procedure. The Meditech MI 111 Universal Surgical Table features can be opted by its mechanism i.e. Lateral tilt, Trendelenburg/ Reverse Trendelenburg, Chair Position, Flex – reflex and Adjustable height functions. The Meditech MI 111 Universal Surgical Table is constructed of stainless steel and other high quality materials. The table is equipped with a large translucent top, for C-arm use.



### Features

- Specialized for orthopaedic, neurology, ERCP, laparoscopy and urology.
  - The MI 111 provides excellent C-Arm and surgical access.
  - Slimline column offset relative to tabletop and base.
  - Furnished with hard, robust stainless steel shell which provides the sense of dirt free & sanitized operations.
  - Smooth height adjustment of the table can be achieved by using the foot pedal.
  - Trendelenburg, lateral tilt, flex-reflex, chair operated manually.
- Interchangeable general sections (head/leg).



CAPITAL MEDICAL EQUIPMENTS SALES & SERVICES

**MediworldIndia**

contact@mediworldindia.com  
mediworld21@gmail.com  
+91 98900 67867 or +91 8668334988



CAPITAL MEDICAL EQUIPMENTS SALES & SERVICES

**MediworldIndia**

contact@mediworldindia.com  
mediworld21@gmail.com  
+91 98900 67867 or +91 8668334988

## Technical Specification

<b>Top dimension</b>	L 1905 x W 533 mm
<b>Height adjustment</b>	750 mm – 1000 mm
<b>Trendelenburg / Reverse</b>	30° / 25°
<b>Lateral tilt</b>	20° / 20°
<b>Kidney elevator</b>	150 mm
<b>Back Rest (up / down)</b>	80° / 25°
<b>Leg Rest (up / down)</b>	15° / 90°
<b>Head Rest (up / down)</b>	20° / 60°
<b>Patient Weight Capacity</b>	150 kgs (330 lbs)



CAPITAL MEDICAL EQUIPMENTS SALES & SERVICES

**MediworldIndia**

contact@mediworldindia.com  
mediworld21@gmail.com  
+91 98900 67867 or +91 8668334988