



CAPITAL MEDICAL EQUIPMENTS SALES & SERVICES

MediworldIndia

contact@mediworldindia.com
mediworld21@gmail.com
+91 98900 67867 or +91 8668334988

Electro Manual OT Table with Dual Operations

Features

- Electro Manual drives with noiseless, jerk less and smooth positioning.
- Table positions i.e. **Height adjustment, Trendelenburg, Reverse Trendelenburg, Lateral Tilt, Back Section/Chair Position, Flex-Reflex are controlled by dual operations i.e. Corded Handset & Manual Override with graphical icons.**
- Zero Position to make table neutral by simply pressing '0' positioning button.
- Remote has **LED backlit** with symbolic position figures to use even in darkness.
- **Safety key** to block motor system in critical operations.
- Manually gear operated Kidney Elevator.
- **Table top is divided into 5 sections** providing numerous articulation possibilities.
- Radiolucent table top is **C-ARM & X-RAY compatible** over the entire length of the table.
- No cross bar on the table top making it very suitable for fluoroscopy.
- The base of the table is scientifically designed "T" shaped SS covered heavy base providing optimum foot room for comfortable posture of the surgeons.
- Offset column provide unrestricted access of C-ARM.
- Detachable-Interchangeable head and leg sections provide further versatility.
- Integrated weight compensation by manual adjustment of the head and leg sections.
- **Non-corrosive stainless steel 304G construction including table top frames, accessories rail, base cover and column telescopic covers.**
- Four easy glide, antistatic castors, two of which are swiveling, ensure mobility.
- **50mm decompressed mattress** with greater comfort & easy to clean.
- **Pedal operated braking system** secures the table in its place while in use.



CAPITAL MEDICAL EQUIPMENTS SALES & SERVICES

MediworldIndia

contact@mediworldindia.com
mediworld21@gmail.com
+91 98900 67867 or +91 8668334988



CAPITAL MEDICAL EQUIPMENTS SALES & SERVICES

MediworldIndia

contact@mediworldindia.com
mediworld21@gmail.com
+91 98900 67867 or +91 8668334988

Optional

- **Override Panel/Standby Controller:** Additional control panel on the table column cover for continue operation in case of remote failure.
- **Zero Position** to make table neutral by simply pressing '0' positioning button.
- **Memory Function** – The microprocessor based control system is fully programmable and can store upto two user-programmable preset positions in its memory that can be recalled anytime by simply pressing M1 or M2 button.
- **Remote** has LED backlit with symbolic position figures to use even in darkness.
- **Safety key** to block motor system in critical operations.
- **Individual Locking Function.**

Technical Constraints

Top dimension	L 1980 x W 533 mm
Height adjustment with pad	750 mm – 1050 mm
Trendelenburg / Reverse	30° / 25°
Lateral tilt left/right	20° / 20°
Kidney elevator	150 mm
Back Rest (upward / downward)	80° / 25°
Leg Rest (upward / downward)	15° / 90°
Head Rest (upward / downward)	20° / 90°
Power supply	24VDC
Patient Weight Capacity	200 kgs (440 lbs)



CAPITAL MEDICAL EQUIPMENTS SALES & SERVICES

MediworldIndia

contact@mediworldindia.com
mediworld21@gmail.com
+91 98900 67867 or +91 8668334988