



CAPITAL MEDICAL EQUIPMENTS SALES & SERVICES

MediworldIndia

contact@mediworldindia.com
mediworld21@gmail.com
+91 98900 67867 or +91 8668334988

Electromatic Operation Table

Description

The Meditech MI 101 Universal Surgical Bed is a mobile, electrical operated surgical table designed to support Cardiology, Orthopedic, Neurology, Gynecology, Urology and ERCP Procedures. The table is equipped with a large translucent top providing maximum radiological access & C-ARM functions. The Meditech MI 101 Universal Surgical Table is powered by Electric power system



Features

- Specialized for cardiology, orthopaedic, neurology, ERCP, laparoscopy and urology.
- The MI 101 offers outstanding imaging access and functionality, with a full length translucent top.
- Offset positioning allows perfect access of C-Arm and clear access for surgeon.
- Interchangeable head and leg sections offer optimum patient positioning.
- The table has been made from acid proof, antirusted stainless steel 304G (the face of the table is easy to clean and immune to disinfection agents).
- Corded handset for simple, intuitive control (height, trendelenburg, lateral tilt, flex-reflex, chair).

Optional

- **Override Panel/Standby Controller:** Additional control panel on the table column cover for continue operation in case of remote failure.
- **Zero Position** to make table neutral by simply pressing '0' positioning button.
- **Memory Function** – The microprocessor based control system is fully programmable and can store upto **two user-programmable preset positions** in its memory that can be recalled anytime by simply pressing M1 or M2 button.
- Remote has **LED backlit** with symbolic position figures to use even in darkness.
- **Safety key** to block motor system in critical operations.
- Individual Locking Function.



CAPITAL MEDICAL EQUIPMENTS SALES & SERVICES

MediworldIndia

contact@mediworldindia.com
mediworld21@gmail.com
+91 98900 67867 or +91 8668334988



MediworldIndia

contact@mediworldindia.com
mediworld21@gmail.com
+91 98900 67867 or +91 8668334988

CAPITAL MEDICAL EQUIPMENTS SALES & SERVICES

Technical Specification

Top dimension	L 1930 x W 533 mm
Height adjustment	750 mm – 1000 mm
Trendelenburg / Reverse	30° / 25°
Lateral tilt	20° / 20°
Kidney elevator	150 mm
Back Rest (up / down)	80° / 25°
Leg Rest (up / down)	15° / 90°
Head Rest (up / down)	20° / 60°
Power supply	24VDC
Battery Backup	3 hours (optional)
Patient Weight Capacity	200 kgs (440 lbs)



MediworldIndia

contact@mediworldindia.com
mediworld21@gmail.com
+91 98900 67867 or +91 8668334988

CAPITAL MEDICAL EQUIPMENTS SALES & SERVICES