



CAPITAL MEDICAL EQUIPMENTS SALES & SERVICES

MediworldIndia

contact@mediworldindia.com
mediworld21@gmail.com
+91 98900 67867 or +91 8668334988

Horizontal Cylindrical Steam Sterilizer

Salient Feature

End Ring Hydraulically formed with out joint to sustain Required pressure

- All other joints Argon Arc welded.
- Pressure lock device.
- Pressure lock device for the door.
- Chamber, Jacket, Door Locking Ring & Back Plate is made from Stainless steel.
- Magnetic Float Switch to Protect Water Immersion Heater from low Water Level.
- Two spring Relief Safety valve to Prevent Un controlled Pressure.
- All Measuring Components (Pressure gauge & Temperature Gauge) are calibrated.
- Vacuum system incorporated.
- fibre Glass wool Insulation With stainless Steel Outer Cover.
- Fulfill all Requirements of 3829 Part 1 Part 2.
- Power 9 KW. 440 V,3 – Phase & Neutral.
- First Cycle 90 Minutes.
- Subsequent 40 minutes.
- Best suited for Rapid Sterilizations of Instruments for Crdiac, Ortho, Neuro & Urology
- Boiler / Pipe steam any Combination for all.
- PC /PLC based Auto Control Panel. (optional)
- Vacuum pump for drying.(optional)



CAPITAL MEDICAL EQUIPMENTS SALES & SERVICES

MediworldIndia

contact@mediworldindia.com
mediworld21@gmail.com
+91 98900 67867 or +91 8668334988