



CAPITAL MEDICAL EQUIPMENTS SALES & SERVICES

MediworldIndia

contact@mediworldindia.com
mediworld21@gmail.com
+91 98900 67867 or +91 8668334988

Horizontal Rectangular Steam Sterilizer

Description

The **Horizontal Cylindrical Steam Sterilizer** are made of heavy gauge stainless steel with leak proof argon arc welding. The radial locking lid is made of thick S.S. plate.



Salient Features

- Best suited for parenteral solutions, injectables, mattresses And hospital dressings.
- Unipoint control station for all phases of the sterilizing cycle In smaller size.
- Operating Temperature 121°C or higher 134 °C.
- All necessary safety precautions namely Pressure locking Door, Spring loaded Safety valve, dead weight Safety valve, Etc.. Incorporated.
- Loading carriage and trolley as standard accessories for easy transfer of load.
- Choice of Single Door, Double Door, Recessed or Panel Models.
- Various material combination like Mild Steel, Stainless Steel 304/316 is available.
- Steam Trap is incorporated in discharge line for fast and perfect elimination of air.

Optional Features

- Single Door / Double Door.
- Flush Mounted.
- Digital Temperature Recorder.
- Rapid Heating
- Rapid Cooling (Air & Water). High Vacuum Pump.
- High Vacuum Pump
- PC(Computer)/ PLC based Auto Control Panel.
- Air Oven Steam Pulsing Pressure System



CAPITAL MEDICAL EQUIPMENTS SALES & SERVICES

MediworldIndia

contact@mediworldindia.com
mediworld21@gmail.com
+91 98900 67867 or +91 8668334988