



CAPITAL MEDICAL EQUIPMENTS SALES & SERVICES

MediworldIndia

contact@mediworldindia.com
mediworld21@gmail.com
+91 98900 67867 or +91 8668334988

Jupiter Anaesthesia Machine

Features

- Deep and spacious drawer units & convenient footrest
- Wide monitoring shelf at eye level with securing provision
- Convenient swivel common gas outlet
- Hook for breathing circuit
- Wipe-clean stainless steel table top
- Auxiliary high pressure Oxygen outlets.
- Extended platform at rear for two 10 Ltr water capacity cylinders* (optional)
- Open / Close circuit switch control
- Pin indexed yokes for O₂ & N₂O
- Convenient mounting of circle absorber
- Choice of backbar for SELECTATEC or Cagemount vaporisers
- Active / passive anaesthetic gas scavenging interface and suction (optional)
- Integrated rear bumper/ side handle for improved maneuverability
- Large, roller bearing, anti-static castors with front wheel brakes



CAPITAL MEDICAL EQUIPMENTS SALES & SERVICES

MediworldIndia

contact@mediworldindia.com
mediworld21@gmail.com
+91 98900 67867 or +91 8668334988