

# Jupiter Plus Anaesthesia Workstation

## Features

- Rust free framework
- Deep spacious drawer units
- Dual cascaded rotameters for O<sub>2</sub> and N<sub>2</sub>O as standard
- Wide monitoring shelf at eye level
- Convenient common gas outlet
- Hook for breathing circuit
- Open-Close circuit switch control
- Convenient mounting of circle absorber
- Twin Selectatec mountings for 2 Agent Specific Vaporisers
- Pin Indexed yokes Two for O<sub>2</sub> and One for N<sub>2</sub>O
- Anaesthetic gas scavenging port
- Large, roller bearing, anti-static castors with front wheel brakes
- Electrical outlets for connecting monitors
- Central pipeline connections with gas specific NIST connectors

