



CAPITAL MEDICAL EQUIPMENTS SALES & SERVICES

**MediworldIndia**

contact@mediworldindia.com  
mediworld21@gmail.com  
+91 98900 67867 or +91 8668334988

## Laparo Surgery OT Table

### Description

Ergonomically design, Laparo surgery OT Table is manufactured using high grade materials. It is engrossed with automatic, remote controlled C-arm imaging equipment. Furthermore, it readily meets the quality standards and requirements of our customers. The table top is large is size to allow maximum radiological access and C-arm functions.



### Technical Specification

<b>Brand</b>	Meditech
<b>Color</b>	White
<b>Kidney Elevator</b>	150 mm
<b>Top Dimension</b>	L 1905 x W 533 mm
<b>Lateral tilt</b>	20 / 20
<b>Wheel</b>	Yes
<b>Back Rest</b>	80 / 25(up / down)



CAPITAL MEDICAL EQUIPMENTS SALES & SERVICES

**MediworldIndia**

contact@mediworldindia.com  
mediworld21@gmail.com  
+91 98900 67867 or +91 8668334988