



CAPITAL MEDICAL EQUIPMENTS SALES & SERVICES

**MediworldIndia**

contact@mediworldindia.com  
mediworld21@gmail.com  
+91 98900 67867 or +91 8668334988

## Remote Operated Obstetric OT Table

### Features

- Electromatic drives with noiseless, jerk less and smooth positioning.
- Three section upholstered table top ensures maximum comfort to the patient.
- Table positions i.e. Height adjustment, Trendelenburg, Reverse Trendelenburg, Back Section/Chair Position are controlled by corded handset with graphical icons.
- Detachable leg sections offer perineal access & allow to attach drainage tray.
- Integrated weight compensation by manual adjustment of leg section.
- Large size stainless steel drainage tray
- Non-corrosive stainless steel 304G construction including table top frames, accessories rail, base cover and column telescopic covers.
- Easy for cleaning & hygiene.
- The powered movements are activated by Imported Linear Actuators.
- Four easy glide, antistatic castors, two of which are swivelling, ensure mobility.
- Pedal operated central braking system secures the table in its place while in use.

### Technical Constraints

<b>Top dimension</b>	L 1750 x W 600 mm
<b>Height adjustment with pad</b>	700 mm – 950 mm
<b>Trendelenburg / Reverse Trend</b>	15° / 15°
<b>Back Rest (upward / downward)</b>	60° / 0°
<b>Leg Rest (upward / downward)</b>	15° / 90°
<b>Power supply</b>	24VDC
<b>Patient Weight Capacity</b>	150 kgs (330 lbs)



CAPITAL MEDICAL EQUIPMENTS SALES & SERVICES

**MediworldIndia**

contact@mediworldindia.com  
mediworld21@gmail.com  
+91 98900 67867 or +91 8668334988