



CAPITAL MEDICAL EQUIPMENTS SALES & SERVICES

MediworldIndia

contact@mediworldindia.com
mediworld21@gmail.com
+91 98900 67867 or +91 8668334988

Three Pin Head Mayfield Fixrator

Description

We offer a supreme range of Three Pin Head Mayfield Fixrator, which is suitable for brain surgery. This Neuro attachment has been developed & manufactured with innovative quality to meet the demanding needs of Neurosurgeons



Features

- The skull clamp has the additional feature of a swiveling C-Arm.
- This swiveling capacity greatly facilitates positioning the pins on the head around nerves, vessels, and thin bone in the temporal region.
- For better result, the arm is locked by control knob and the width is adjusted by a ratchet mechanism.
- The three cone-shaped pins available in adult and child sizes are interchangeable and may be removed for sterilization.



CAPITAL MEDICAL EQUIPMENTS SALES & SERVICES

MediworldIndia

contact@mediworldindia.com
mediworld21@gmail.com
+91 98900 67867 or +91 8668334988